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Uro-/ Urogynäkologischer Studien

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1. Plättchen-reiches Plasma: Grundlagen	3
1.1. Standardisiertes PRP: Gomri et. al. (2022).....	6
2. Studien zu PRP im Uro-/ Urogynäkologischen Kontext	19
2.1. Urologische Studie: Zaghloul A. et. al. (2022)	19
2.2. Urologische Studie: Poullos, E. et. al. (2023)	20
2.3. Urogynäkologische Studie: Prodromidou, A. et. al. (2021)	22
2.4. Urologische Studie: Yáñez-Castillo et. al. (2022)	24
2.5. Urologische Studie: Zugail, A. S. et. al. (2023)	25
2.6. Gynäkologische Studie: Saleh, M. D. & Abdelghani R. (2022)	26
2.7. Gynäkologische Studie: Cheng-Yu et. al. (2016).....	27
2.8. Gynäkologische Studie: Castellani et. al. (2021)	29
2.9. Gynäkologische Studie: Tognazzo et. al. (2023)	30
2.10. Gynäkologische Studie: Aretaieion Hospital et. al. (2019)	32
3. Weitere relevante Erkenntnisse.....	34
3.1. Virag R. et. al. (2018)	34
3.2. Sexualmedizin	36

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1.1. Standardisiertes PRP: Gomri et. al. (2022)

Open Access Review

Swiss Medical Devices for Autologous Regenerative Medicine: From Innovation to Clinical Validation

by  Farid Gomri ^{1,†}  Solange Vischer ^{1,†}  Antoine Turzi ¹ and  Sarah Berndt ^{1,2,*} 

¹ Regen Lab SA, 1052 Le Mont-sur-Lausanne, Switzerland

² Department of Plastic, Reconstructive and Aesthetic Surgery, Faculty of Medicine, Geneva University Hospitals, 1205 Geneva, Switzerland

* Author to whom correspondence should be addressed.

† These authors contributed equally to this work.

Pharmaceutics **2022**, *14*(8), 1617; <https://doi.org/10.3390/pharmaceutics14081617>

Received: 11 April 2022 / Revised: 15 July 2022 / Accepted: 29 July 2022 / Published: 2 August 2022

(This article belongs to the Special Issue *Advanced Pharmaceutical Science and Technology in Switzerland*)

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Table 2. Published clinical studies in which biologics manufactured with Regen Lab devices were successfully used.

Therapeutic Domain	Pathologies	Studies	Treated Patients
Sports medicine	Osteoarthritis and other cartilage pathologies	37	2053
	Tendinopathies	15	506
	Muscle	3	18
	Carpal tunnel	6	215
Orthopedic surgery	Bone reconstruction, cartilage and tendon surgery	16	649
	Spine surgery	3	46
Dermatology	Skin care	14	353
	Alopecia	10	240
	Chronic wounds	12	238
	Plastic and reconstructive surgery	7	178
	Cardiovascular surgery wounds	3	425
	Urology	4	122
	Gynecology	8	105
	Penile dysfunction	3	151
	Fertility	3	67
	Ophthalmology	10	170
	Otominolaryngology	7	103
	Dental/maxillofacial surgery	7	200
	Total	168	5929

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2.1. Urologische Studie: Zaghloul A. et. al. (2022)

Received: 21 December 2021 | Revised: 8 February 2022 | Accepted: 6 March 2022
DOI: 10.1111/and.14421

ORIGINAL ARTICLE

ANDROLOGIA WILEY

Assessment of the intracavernosal injection platelet-rich plasma in addition to daily oral tadalafil intake in diabetic patients with erectile dysfunction non-responding to on-demand oral PDE5 inhibitors

Ahmed S. Zaghloul | Abdelrahman M. El-Nashaar | Shady Z. Said | Ihab A. Osman | Taymour Mostafa

Department of Andrology, Sexology & STIs, Faculty of Medicine, Cairo University, Cairo, Egypt

Correspondence
Ahmed S. Zaghloul, Department of Andrology, Sexology & STIs, Faculty of Medicine, Cairo University, Cairo 11562, Egypt.
Email: ahmadzag@gmail.com

Abstract

This study aimed to evaluate the effectiveness of ICI of platelet-rich plasma (PRP) in addition to daily oral tadalafil intake in diabetic erectile dysfunction (ED) patients non-responding to PDE5 inhibitors. Overall, 48 patients complaining of ED non-responding to on-demand PDE5 inhibitors were allocated into 2 equal groups, diabetics and non-diabetics that were given a daily dose of 5 mg tadalafil plus vardenafil 20 mg on demand during the study besides being subjected to 3 doses of ICI of PRP, 4 weeks apart. Responses to on-demand PDE5 inhibitors, International index of erectile function-5 (IIEF-5) score, erection hardness scores (EHS) and pharmacodynamic duplex studies were assessed. After PRP injections, 33% and 50% of cases were satisfied with on-demand PDE5 inhibitors, respectively, whereas 41% and 66% of them showed improved EHS response. Compared with baseline scores, the mean IIEF-5 scores were significantly improved after PRP therapy in the diabetic ED group (12.1 vs. 8.04, $p = 0.003$) as well as in the non-diabetic ED group (14.8 vs. 10.2, $p = 0.001$) linked to pharmaco-penile duplex readings. Both good and fair diabetic control exhibited significant responses to ICI therapy of PRP compared with bad controlled cases. The significant improvement included; the IIEF-5 score increase (86.7%, 126% vs. 16.1%), improved EHS as well as penile duplex readings. Baseline HbA1C demonstrated a significant negative correlation with IIEF-5 score before ($p = 0.019$) and after PRP therapy ($p = 0.002$) respectively. It could be concluded that ICI of PRP could be an effective therapy for treating ED patients non-responding to on-demand oral PDE5 treatment.

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2.2. Urologische Studie: Poullos, E. et. al. (2023)

> Sex Med Rev. 2023 Sep 27;11(4):359-368. doi: 10.1093/sxmrev/qead027.

Platelet-rich plasma for the treatment of erectile dysfunction: a systematic review of preclinical and clinical studies

Evangelos Poullos¹, Ioannis Mykoniatis^{1,2}, Nikolaos Pyrgidis^{2,3}, Dimitrios Kalyvianakis^{1,2}, Dimitrios Hatzichristou^{1,2}

Affiliations + expand

PMID: 37528499 DOI: 10.1093/sxmrev/qead027

Abstract

Introduction: There has been tremendous growth in regenerative medicine during the last decade. For erectile dysfunction (ED), after the inclusion of low-intensity shockwave therapy as a treatment modality for ED management by the European Association of Urology sexual health guidelines, intracavernosal injection of platelet-rich plasma (PRP) has gained popularity between urologists and patients as a novel ED therapeutic modality with initial promising results. However, limited clinical data exist regarding efficacy and safety in patients with ED. Furthermore, despite numerous preclinical studies in other tissues and organs, the mechanism of action for restoring erectile function remains undetermined.

Objectives: This systematic review aims to present the current status of preclinical and clinical evidence regarding the use of PRP as treatment option for ED.

Methods: A systematic literature search was conducted using PubMed, Cochrane, and ScienceDirect databases, until February 2023 for studies exploring the effect of PRP on ED.

Results: We identified 517 articles, 23 of which were included in this review. These were 7 preclinical (of which 1 was a comparative trial and 6 were placebo-controlled randomized controlled trials) and 16 clinical studies (of which 1 was a comparative trial, 5 were randomized trials, and 2 were placebo-controlled randomized controlled trials). Preclinical data support the regenerative role of PRP in erectile tissue, in accordance with existing evidence in other tissues. Randomized clinical studies, as well as the first 2 available randomized, placebo-controlled clinical trials, showed promising efficacy and a lack of any adverse events.

Conclusion: As PRP for ED is widely used worldwide, there is an urgent need for high-quality studies with long-term follow-up. Standardization of research protocols, especially on the quality of PRP preparation, is also needed.

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2.3. Urogynäkologische Studie: Prodromidou, A. et. al. (2021)

The Emerging Role on the Use of Platelet-Rich Plasma Products in the Management of Urogynaecological Disorders

Surgical Innovation
2021, Vol. 0(0) 1–8
© The Author(s) 2021
Article reuse guidelines:
sagepub.com/journals-permissions
DOI: 10.1177/15533506211014848
journals.sagepub.com/home/sri

Anastasia Prodromidou, MD, PhD¹ , Dimitrios Zacharakis, MD, PhD¹, Stavros Athanasiou, MD, PhD¹, Athanasios Protopapas, MD, PhD¹, Lina Michala, MD, PhD¹, Nikolaos Kathopoulos, MD, PhD¹ , and Themis Grigoriadis, MD, PhD¹

Abstract

Background: The regenerative efficacy of platelet-derived products has been recently investigated in the treatment of pelvic floor disorders (PFDs). We aimed to synthesize the current evidence of platelet-rich plasma (PRP) products used in urogynaecological disorders including vaginal atrophy, pelvic organ prolapse (POP), urinary incontinence, vaginal fistulas and vaginal mesh exposure. **Methods:** A meticulous search of the currently available literature on the use of PRP for the management of PFDs was performed using 3 electronic databases. **Results:** PRP could be a feasible alternative modality for the management of vaginal atrophy with favourable outcomes in vaginal atrophy parameters and patients' satisfaction, especially when hormone therapy is contraindicated. In patients with POP, an increase in collagen concentration after PRP application was observed while the use of PRP resulted in improvement of stress urinary incontinence symptoms. A considerable proportion of vesicovaginal fistulas were treated after application of PRP-based injections. **Conclusions:** There is only limited evidence of the use of PRP for PFDs. Platelet-rich plasma appears to be a promising, easy to apply, cost-effective and feasible alternative therapeutic modality for the management of various urogynaecological disorders. Future randomized trials are needed to confirm the efficacy of PRP in the treatment of urogynaecological disorders.

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Table 1. Characteristics of the included studies for the use of PRP in vaginal atrophy.

Year; author	Type of study	Country	Type of gynaecological disorder/procedure	Inclusion criteria	Assessed parameters	PRP preparation technique	Dose of PRP injection	Injected areas
2018; Hersant	Phase II pilot PS	France	Vulvo-vaginal laxity/PRP-HA injections	>18 years; history of hormone-dependent breast cancer in remission >5 years after surgery; VHI<15; no vulvo-vaginal inflammation or infection; no history of recurrent vaginal herpes; no lichen sclerosis; no history or genital cancer	VHI score>15; vaginal pH; efficacy of treatment; sexual quality; pain (VAS score)	RegenKit BCT-HA (a sterile, nonpyrogenic tube designed for use in preparing a mixture of PRP and HA (40 mg [2% w/v] of non-cross-linked HA per tube, 1550 KDa)); centrifugation 1500 g for 5 minutes	4 mL of PRP-HA (2 mL PRP mixed with 2 mL HA) (.1 mL/point) Injections every 5 mm with 27G calibre needle and a 1-mL syringe	Into the vestibule and the 1st 3 cm of the posterior vagina and 2 cm posterior wall of the introitus using a speculum or by laterally opening the vaginal walls with the fingers
2017; Kim	CR	Korea	Vaginal atrophy and lichen sclerosis/filling with autologous fat and PRP	N/A	Postoperative photographs	Double spin centrifugation with SmartPrep APC-30 4cc of PRP derived from 30cc of autologous PRP	40cc autologous fat mixed with PRP into 1cc syringes	Subcutaneous layer of the labia majora aseptically via 4 ports
2016; Aguilar	CR	France	Vaginal laxity-sexual dysfunction/ vaginoplasty + PRP with HA	N/A	Sabbatsberg sexual self-rating scale	RegenKit BCT-HA 3 sterile tubes (4 mL/tube). Each contained a gel of HA, and inert cell-selector gel, a liquid anticoagulant and allowed the preparation of 4 mL of PRP-HA mixture; centrifugation: 1500 g for 5 minutes (cellular elements settle on the surface) manual homogenization PRP and HA; harvesting fat with Coleman's technique centrifuged at 1500 g for 1 minute	4 mL of PRP-HA (2 mL PRP mixed with 2 mL HA)	16 mL of fat cells in the posterior vaginal wall (2 linear anteroposterior tracts) and 10 mL PRP-HA in the perineal raphe (episiotomy scar) in the vestibular fossa and in the labium minus and majus

Abbreviations: CS = case report; HA = hyaluronic acid; N/A = not available; PS = prospective; PRP = platelet-rich plasma; VAS = visual analogue scale; VHI = vaginal health index.

Table 2. Characteristics of the included studies for the use of PRP in POP.

Year; author	Type of study	Country	Type of gynaecologic disorder/treatment	Inclusion criteria	Assessed parameters	PRP preparation	Dose of PRP	Type of application
2015; Medel	In vitro	Canada	POP/12 weeks after POP repair samples from 6 patients collected to obtain PRP	N/A	Attachment of POP human vaginal fibroblast to 2 meshes	Regen ACR-C kit 8 mL whole blood collected into EDTA-coated tube centrifuge for 5 minutes at 1500 g; collection of the	Coat in 200 µL PRP	Coating of 5 × 5mm squares of absorbable polyglactin mesh (5 samples) and nonabsorbable polypropylene mesh

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Applications of Platelet Rich Plasma in Urology

Yaiza María Yáñez-Castillo ¹, Miguel Ángel Arrabal-Polo ¹

Affiliations + expand

PMID: 36330568 DOI: 10.56434/j.arch.esp.urol.20227508.98

Free article

Abstract

Introduction: Platelet rich plasma is a product obtained from vein whole blood sample, with a high concentration of growth factors, which have been reported to be beneficial for tissue regeneration. In the last few years, several assays that pretend to demonstrate the efficacy of platelet rich plasma in many disciplines, including urology, have been published.

Aim and methods: The aim of this assay is to develop a bibliographic review of the publications available about platelet rich plasma and urology. The search was based on PubMed database.

Conclusions: Platelet rich plasma has demonstrated to be a safe and easy to obtain product, and it has also shown promising results in terms of efficacy in urology. However, it is necessary to conduct large, placebo-controlled, studies that can confirm or deny the effectiveness of platelet rich plasma before it can be used on urological conditions.

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ARTICLE

 Check for updates

Safety and feasibility of percutaneous needle tunneling with platelet-rich plasma injections for Peyronie's disease in the outpatient setting: a pilot study

Ahmed S. Zugail ^{1,2}, Muaath Alshuaibi ^{1,3}, Sandrine Lombion⁴ and Sébastien Beley¹

© The Author(s), under exclusive licence to Springer Nature Limited 2023

The objective of this study is to evaluate the safety and feasibility of the combined simultaneous percutaneous needle tunneling coupled with injection of platelet-rich plasma in the outpatient department for the treatment of Peyronie's disease. This prospective, non-randomized, cohort and preliminary study included patients who underwent this procedure from November 2020 to July 2022. The main outcome was an improvement in penile curvature. Fifty-four patients were enrolled and underwent 6 sessions under local anesthesia followed by vacuum therapy for the treatment of Peyronie's disease in our outpatient unit. The amendment of the curvature angle was significant with a median correction percentage of -44.40% interquartile range (-66.70 to -39.70), [p -value = 0.001 , 95% CI (-29.76 to -18.02)], paired Student's t -test]. The median pre-treatment curvature angle was 45° (40 – 75), and the median post-treatment was 30° (20 – 40). The median score for pain during the procedure was 3 (0 – 4.25) according to a 10-point visual analogic scale. After two hours, 20.37% of patients still had pain but none required any pain medication. 50% of patients had a minor hematoma and 75.93% patients reported penile ecchymosis. A single patient reported an injection site skin infection. In our experience percutaneous needle tunneling with platelet-rich plasma injections for Peyronie's disease in the outpatient setting is a safe, effective, and feasible treatment of penile deformity for PD.

https://www.nature.com/articles/s41443-023-00744-y.epdf?sharing_token=V_29iVjPtSxScWnY46uwf9RgN0jAjWel9jnR3ZoTv0MeKMoRJuSUM7BJpMDr7WzNavl2qQweTJC8mC2n7Z5O4gOwwXBpQguAPsmdkwXpftzHann_IK7nk41euti_ekUzx6Flx_mbYDBxD5c1LrQZ-kxJL2mKB7tnIRQsVKc0LrZw%3D

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2.6. Gynäkologische Studie: Saleh, M. D. & Abdelghani R. (2022)

Received: 20 January 2022 | Accepted: 14 February 2022
DOI: 10.1111/jocd.14873

ORIGINAL ARTICLE

Clinical evaluation of autologous platelet rich plasma injection in postmenopausal vulvovaginal atrophy: A pilot study

Doaa M. Saleh MD^{1,2} | Rania Abdelghani MD^{2,3}

¹Obstetrics and Gynecology, Faculty of Medicine for Girls, Al-Azhar University, Cairo, Egypt

²Armed Forces College of Medicine (AFCM), Cairo, Egypt

³Dermatology and Venereology, Faculty of Medicine for Girls, Al-Azhar University, Cairo, Egypt

Correspondence
Rania Abdelghani, Dermatology and Venereology, Faculty of Medicine for Girls, Al-Azhar University, 3 Al-Fedaa street, Nasr City, 11371 Cairo, Egypt.
Emails: rania_abdelghani@hotmail.com; ranyaabdelghani.med@azhar.edu.eg; rania_abdelghani@fcm.edu.eg

Funding information
No funding.

Abstract

Background: There is a lack of published data investigating injection of autologous platelet rich plasma (A-PRP) alone in the treatment of postmenopausal VVA.

Objectives: In this pilot study, we aimed to investigate the safety and efficacy of injection of A-PRP alone in postmenopausal VVA in women without the history of cancer breast to explore its utility as a hormone free therapy for postmenopausal VVA and for vulvovaginal rejuvenation.

Methods: In this pilot study, 47 women with postmenopausal VVA were included. Vulvovaginal condition was evaluated at the baseline by vaginal health index (VHI). Impact of VVA on quality of life and sexual life was evaluated at the baseline by vulvovaginal symptom questionnaire (VSQ). Treatment protocol was 2 sessions of A-PRP injection with 1 month interval. Response was evaluated 1 month after the last session by VHI and VSQ. Side effects were also evaluated.

Results: Postmenopausal VVA was significantly improved by A-PRP injection as indicated by significant improvement of total VHI score and its items at 1 month post-treatment (p value <0.001). Moreover, there was significant improvement of burning, hurting, being irritated, being dry, discharge, desire to be intimate, sexual relationships, pain during sexual activity, and dryness during sexual activity at 1 month post-treatment as indicated by VSQ (p value = 0.045 for being dry and <0.001 for other items).

Conclusions: Autologous platelet rich plasma injection is safe and effective as minimally invasive monotherapy for postmenopausal VVA without history of cancer breast and hence for vulvovaginal rejuvenation.

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2.7. Gynäkologische Studie: Cheng-Yu et. al. (2016)

A pilot study: effectiveness of local injection of autologous platelet-rich plasma in treating women with stress urinary incontinence

Cheng-Yu Long ^{# 1 2 3}, Kun-Ling Lin ^{# 2 4}, Chin-Ru Shen ², Chin-Ru Ker ², Yi-Yin Liu ^{1 3}, Zi-Xi Loo ², Hui-Hua Hsiao ⁵, Yung-Chin Lee ^{6 7 8 9}

Affiliations + expand

PMID: 33452330 PMCID: PMC7810841 DOI: 10.1038/s41598-020-80598-2

[Free PMC article](#)

Abstract

The study aims to evaluate the effectiveness of local injection of autologous platelet rich plasma (A-PRP) as a treatment for women suffering from stress urinary incontinence (SUI). In a prospective intervention study, twenty consecutive women suffering from SUI were treated with A-PRP injection at anterior vaginal wall where mid-urethra locates. Self-reported questionnaires were used to measure pre-treatment, 1 month and 6 months post-treatment symptom severity. Secondary outcomes of sexual function and treatment effect sorted by age were analyzed with valid statistical methods. A-PRP is effective in relieving SUI symptoms at both 1 month and 6 months post-treatment without significant adverse reactions reported. It seems to have a trend that treatment success rate with cured and improved symptoms was slightly higher in the younger group, although it did not reach statistical significance ($P = 0.07$). No significant changes in sexual function before and after the treatment were reported by the patients. This pilot study is the first to report A-PRP treatment effect for SUI in women. The result suggested that A-PRP is a considerable treatment option for mild to moderate SUI cases. It also opens up further research opportunities for A-PRP's clinical applications.

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An innovative approach to treating vaginal mesh exposure after abdominal sacral colpopexy: endoscopic resection of mesh and platelet-rich plasma; initial experience in three women

Daniele Castellani¹, Alessandra Valloni², Angela Piccirilli², Giuseppe Paradiso Galatioto², Carlo Vicentini²

Affiliations + expand

PMID: 27655187 DOI: 10.1007/s00192-016-3154-x

Abstract

Introduction and hypothesis: Polypropylene mesh exposure is uncommon after abdominal sacral colpopexy (ASC), but in case of symptomatic vaginal mesh exposure, surgery is needed. When treating it, care must be taken to completely remove the exposed mesh (EM), saving as much vaginal tissue as possible to avoid a subsequent shortened and narrowed vagina. In this video, we present a minimally invasive technique for treating EM after ASC using endoscopic mesh resection and autologous platelet-rich plasma (PRP) technology.

Methods: Three women were referred to our outpatient clinic for vaginal vault mesh exposure after laparoscopic ASC with concomitant hysterectomy. All women underwent endoscopic bipolar PlasmaKinetic resection (BPR) of EM, and PRP gel was delivered in the surgical site to cover the gap left by the resection.

Results: Mean operative time was 39.6 min. Surgery was uneventful in all cases. All women recovered sexual function, and nobody experienced relapsed pelvic organ prolapse at 1-year follow-up.

Conclusions: Our preliminary results show that BPR and PRP are safe, effective, and feasible for treating vaginal mesh exposure with conservation of anatomy and sexual function.

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2.9. Gynäkologische Studie: Tognazzo et. al. (2023)

JOURNAL ARTICLE

Autologous Platelet-Rich Plasma in Clitoral Reconstructive Surgery After Female Genital Mutilation/Cutting: A Pilot Case Study [Get access](#)

Enrico Tognazzo, MD, Sarah Berndt, PhD, Jasmine Abdulcadir, MD, PD, FECSM

Aesthetic Surgery Journal, Volume 43, Issue 3, March 2023, Pages 340-350,
<https://doi.org/10.1093/asj/sjac265>

Published: 18 October 2022 [Article history](#)

“ Cite Permissions Share

Abstract

Background

Postoperative healing after clitoral reconstruction (CR) for female genital mutilation/cutting can be long and painful due to prolonged clitoral re-epithelialization time (up to 3 months). Autologous platelet-rich plasma (A-PRP) might reduce postoperative clitoral epithelialization time and pain.

Objectives

The authors assessed postoperative clitoral re-epithelialization time and pain after intraoperative clitoral administration of A-PRP.

Methods

Five consecutive women underwent CR (Fold's technique) followed by the administration of A-PRP Regen Lab SA (Le Mont-sur-Lausanne, Switzerland) plasma and glue, injected inside and applied above the re-exposed clitoris, respectively. We recorded surgical complications, postoperative clitoral pain (visual analogue scale), painkiller intake, time to complete re-epithelialization, and the experienced subjective changes in sexual response and perception of their own body image referred by women.

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2.10. Gynäkologische Studie: Aretaieion Hospital et. al. (2019)

HJOG 2019, 18 (3), 55-65

The use of Platelet-Rich Plasma in the Gynaecological Clinical Setting. A review.

Vasileios Sanoulis¹, Nikolaos Nikolettos², Nikolaos Vlahos³

¹Aretaieion Hospital, Athens, Greece, ²Demokriteion University of Thrace, Alexandroupolis, Greece, ³Kapodistrian University of Athens, Athens, Greece

Correspondence

Correspondence should be addressed to Vasileios Sanoulis; vasanoulis@yahoo.gr

Abstract

Reports on clinical use of platelet-rich plasma (PRP) have dramatically increased in the last decade. PRP is a natural product where a high level of platelets is concentrated with growth factors concentration 3 to 5 times greater than plasma. In their alpha granules, platelets contain a cocktail of chemical mediators that trigger tissue regeneration, both through the promotion of angiogenesis and the initiation of cell regeneration. Such bioactive factors include platelet derived growth factor (PDGF), transforming growth factor-beta (TGF- β), vascular endothelial growth factor (VEGF), epidermal growth factor (EGF), fibroblast growth factor (FGF), insulin growth factor (IGF). In the gynaecological setting, PRP has been shown to facilitate wound healing and has been proposed as an effective treatment for the treatment of disorders such as Lichen Sclerosus (LS) or cervical ectopy. Evidence of PRP effectiveness in conditions of pelvic organ prolapse and urinary incontinence has also been reported. Recent publications have shown its positive effects in promoting endometrial and follicular growth and gestation in assisted reproduction cycles. PRP administration in the ovarian tissue has resulted in rejuvenation as evaluated by the relevant biochemical markers as well as reactivation of folliculogenesis. Other researchers have showed effective use in improving pregnancy outcomes in patients with repeated implantation failure. For managing repeated failures, several methods have been attempted, but recently PRP appears as effective in inducing endometrial development in patients with thin, poor endometrium. Another possible mechanism of PRP to endometrium is through its anti-inflammatory actions. PRP is an innovative therapeutic modality, as it is simple, easily performed, and effective, with minimal risk as it is an autologous product. The employment of regenerative medicine approaches will continue to grow in the gynaecological setting and the future trials should be carefully designed in order to facilitate the clinical translation.

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Sexual Health Issues

Research Article

Evaluation of the benefit of using a combination of autologous platelet rich-plasma and hyaluronic acid for the treatment of Peyronie's disease

Ronald Virag^{1*}, Hélène Sussman¹, Sandrine Lambion² and Valérie de Fourmestraux³

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³Regenlab SA, Switzerland

Abstract

Background: Platelet-Rich Plasma (PRP) has been successfully used for different medical conditions, to help regenerate fibrotic or inflammatory tissues. Peyronie's disease (PD) could be a target for such treatment. In this study, PRP in combination with hyaluronic acid (HA) is proposed as an innovative treatment for this condition.

Objectives: To assess the effects of repeated intralesional injections inside the tunica albuginea (TA) with a combination of PRP and HA in PD patients.

Design, setting, and participants: 90 patients, naïve of any prior interventional treatment, were injected with 8 ml of PRP combined with HA directly in the diseased area of TA under ultrasound (US) guidance. Injections were repeated 4 times within 2 months. Additional monthly sessions were done, if necessary.

Outcome measurements and statistical analysis: Penile deformation, TA thickening, presence and size of calcifications, and questionnaires evaluating both PD and sexual function. Global results were classified into satisfactory improvement (I+) unsatisfactory (I-) and unchanged (U).

Results and limitations: Angulation and thickening were both significantly improved after 4 sessions and continuously improved with the number of sessions, to achieve a mean final improvement of 39.65 +/-24.83% (p<0.001) (CI 14.26;18.82) and 26.23 +/-10.29% (p<0.001) (CI 0.99;1.22), respectively. Younger patients achieved better results. 73.3% of the patients belonged to I+ group; 67.8% felt an improvement, 62.2% a decrease of their deformation, 46.7% had easier sexual intercourse and 70% judged the treatment as positive. Limitations included the absence of a control group.

Conclusions: Repeated intralesional injections of PRP-HA improve significantly penile deformation and thickening of the TA. This fully natural treatment may stand as a first line treatment of PD.

Patient summary: We evaluated here an innovative treatment of Peyronie's disease using intralesional injections of Platelet-Rich Plasma mixed with hyaluronic acid. We observed a significant improvement of the deformation in 70% of the cases without any serious adverse events.

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Key Results From Various Publications

7 publications were identified including clinical data on 154 patients. Regen PRP has been used for women for the treatment of vaginal mesh exposure, lichen sclerosus and dyspareunia (female sexual dysfunction) and used to stimulate ovarian regeneration and folliculogenesis activation in peri-menopausal women.

Also, Regen PRP have been used on men for the treatment of Peyronie's disease. Generally, PRP treatment was safe; minimal-moderate pain was reported in few patients and this was limited to the first 24 hours post-treatment. No AE was reported. Performance data show that PRP treatment has a beneficial effect for several gynaecology medical indications; also, PRP was shown to be effective in patients who are not responsive to other therapies such as steroids treatment.

In addition, the combination PRP+HA was shown to be effective for the treatment of Peyronie's disease in men, with no serious adverse event reported. These data contribute to the safety and performance profile of PRP prepared with Regen PRP and support its safety and efficacy for gynaecology medical indications.

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These key results described in publications are summarized in the following table:

Data Reference	Study design & follow-up	Medical indication/ pathology	N° of patients	Treatment protocol	Key results on performance	Safety / adverse events (AE)
<p>“Evaluation of the benefit of using a combination of autologous platelet rich-plasma and hyaluronic acid for the treatment of Peyronie’s disease” Virag R., Sussman H., Lambion S., de Fourmestraux V. <i>Sex Health Issues</i>, 2017 doi: 10.15761/SHI.1000102 Volume 1(1): 1-6</p>	<p>Case series Observational study Follow-up: Up to 10 months</p>	<p>Penile deformations and/or curvatures due to plaques palpated in the tunica albuginea (TA), lasting 20 months in average.</p>	<p>90</p>	<p>All patients underwent 3-steps penile local anaesthesia 4 ml of PRP and 4 ml of a combination of PRP/HA were prepared with a RegenBCT tube and a BCT-HA tube, as per IFU. Both preparations were collected in a syringe and mixed together and injected into the TA’s diseased sites under Duplex Ultrasound control. At the end of the procedure, a slightly compressive dressing was applied. All patients underwent 4 treatment administered with 15-day interval. If needed, treatment was repeated each 3 months.</p>	<p>Deformation and thickening of the TA one month after 4 treatment sessions (n = 84) and at the last follow-up visit (n = 54) showed a significant decrease. The mean final improvement was respectively 16.54 ± 0.51° and 1.11 ± 0.52 mm with a mean gain percentage of $39.65 \pm 24.83\%$ and $26.22 \pm 10.29\%$. Presence of calcifications (n=28) did not change significantly the results but a decrease was detected by DUS examination in 6 patients. Peyronie’s disease questionnaire was significantly improved in all domains. Globally, 62.2% of the patients felt a decrease in their deformation, 43.3% judged that their erections were better, 46.7% thought that sexual activity was easier and</p>	<p>No serious adverse events were recorded. Ecchymosis was present in 16.7% of the patients and marked hematoma in 10%. Two patients suffered from transient hypotension after local anesthesia. One patient stopped the treatment after 4 sessions because of worsening symptoms.</p>
<p>“An innovative approach to treating vaginal mesh exposure after abdominal sacral colpopexy: endoscopic resection of mesh and platelet-rich plasma; initial experience in three women” Castellani D, Valloni A, Piccirilli A, Paradiso Galatioto G, Vicentini C. <i>Int Urogynecol</i> 2017; 28: 325</p>	<p>Case reports Follow-up: 1 year</p>	<p>Vaginal vault mesh exposure after laparoscopic ASC with concomitant hysterectomy (dyspareunia, vaginal bleeding/discharge, perineal pain)</p>	<p>3</p>	<p>Transvaginal endoscopic resection of the extruding prosthetic fragment was conducted by using Iglesias resectoscope and Gyrus device. PRP and autologous thrombin serum were prepared according to IFU. PRP was infiltrated along the margins of the breach created by the resection with a 25G spinal needle. PRP and Autologous Thrombin Serum (ATS) gel was infiltrated along the margins of the breach created by the resection.</p>	<p>67.8% considered that the treatment improved their initial condition. 30 days post-PRP injections, physical urogynecological and vaginocopy examinations showed tissue response to the treatment. Full recovery was obtained already 6-month post injections. All women recovered sexual function, and nobody experienced relapsed pelvic organ prolapse at 1- year follow-up. Bipolar plasmakinetic resection (BPR) combined with PRP is an effective technique for treating vaginal mesh exposure with conservation of anatomical results and sexual function. PRP represents an alternative method to suture stitches, being decisive in solving problems of responsiveness in hysterectomy surgeries</p>	<p>No AE reported. No complication observed during surgery.</p>

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<p>"Use of Platelet-rich Plasma for vulvovaginal autoimmune conditions like Lichen sclerosus" Behnia-Willison F, Pour NR, Mohamadi B <i>et al.</i> <i>Plastic and reconstructive surgery Global open</i> 2016; 4: e1124</p>	<p>Case report</p> <p>Follow-up: 12 months</p>	<p>Lichen sclerosus unresponsive to steroid treatment</p>	<p>28</p>	<p>PRP injected under local anaesthesia (lignocaine, 23%; tetracaine, 7%) to any affected areas of the external genitalia, including the labia majora, labia minora, clitoris, and clitoral hood. The injection was carried out using a 27-gauge needle in a fanning motion to break the scar and fibrotic tissue and retrograde injection of PRP in the tissue.</p> <p>Patients received 3 PRP treatments, 4 to 6 weeks apart and again at 12 months</p>	<p>Based on the Australian Pelvic Floor Questionnaire 24 and colposcopy results, almost patients showed clinical improvement in the size of their lesions, and in 8 cases, lesions totally disappeared after treatment. Symptoms disappeared in 15 of the 28 patients after treatment, with no need for further steroid therapy in 23 patients. 13 women experienced partial symptom relief</p>	<p>Patients reported minimal to moderate pain. During the 24 hours after the procedure, 26 patients (92.9%) reported pain scores of 2 to 3; the remaining patients reported scores of 5 and 7, respectively. No cases of infection, bleeding, hematoma, or other adverse outcomes.</p>
<p>"Ovarian rejuvenation and folliculogenesis reactivation in perimenopausal women after autologous platelet rich plasma treatment" Pantos K, Nitsos N, Kokkali G, <i>et al.</i> <i>In: ESHRE (European Society of Human Reproduction and Embryology) Annual Meeting. Helsinki, Finland; 2016</i></p>	<p>Case reports</p> <p>Follow-up : Not available</p>	<p><i>In vitro</i> ovarian rejuvenation and subsequent folliculogenesis reactivation in perimenopausal women</p>	<p>8</p>	<p>PRP was infused into the ovaries using a transvaginal ultrasound-guided injection technique. When a follicle of >16 mm was observed, ovulation triggering was achieved with 5000 IU of hCG and follicle aspiration was performed 32 hours later by the transvaginal route</p>	<p>The successful ovarian rejuvenation was confirmed by the menstrual cycle restoration 13 months after the ovarian PRP treatment. Oocyte retrieval was successful in all cases, resulting in 2.50±0.71 follicles of 15.20±2.05 mm diameter, 1.50±0.71 oocytes and 1.50±0.71 MII oocytes. All mature oocytes were inseminated by ICSI and the 1.50±0.71 resultant embryos were cryopreserved at 2pn stage until transfer</p>	<p>Not formally assessed.</p>
<p>"A new treatment of lapeyronie's disease by local injections of plasma rich platelets (PRP) and hyaluronic acid. Preliminary results" Ronald Virag, H Sussman, B Lobel <i>www.acad-chirurgie.fr (online version) 1634-0647 - © 2014 Académie nationale de chirurgie</i></p>	<p>Case series</p> <p>Observational study</p>	<p>Lapeyronie's Disease, including 8 patients with curvatures angles >50° and 5 patients with deformities such as</p>	<p>13</p>	<p>8 ml of the PRP+HA combination (prepared with 2 Cellular Matrix BCT-HA tubes, as per IFU) were injected with 18-gauge or 25-gauge needle. All intra-lesional injections were performed under ultrasound-guidance; local anesthesia was applied before treatment administration</p> <p>All patients received 4 treatment sessions administered within 2 months.</p>	<p>The combination of PRP and HA was effective for the majority of patients presenting Peyronie's disease.</p> <p>10 Patients (77%) got improvements with a 30% of decrease in penis curvature. Also, density and size of plaques was decreased in 7 patients (53%).</p> <p>US-images showed five reduction in echogenicity and plaques volume and reduction in density of albuginea images.</p>	<p>Treatment was safe, one patient developed a hematoma, probably related to dorsal vein puncture of which was solved with local ointment.</p>
		<p>strictures and/or deviation <30°</p>				
<p>"A pilot study of the effect of localized injections of autologous platelet rich plasma (prp) for the treatment of female sexual dysfunction" Runels C, Melnick H, Debourbon E, Roy L. <i>J Women's Health Care</i> 2014; 3-4</p>	<p>Case reports</p> <p>Follow-up: 12 weeks</p>	<p>Female sexual dysfunction (dyspareunia)</p>	<p>11</p>	<p>PRP was activated with 0.5 ml of 10% CaCl₂ to form a platelet-rich fibrin matrix (PRFM).</p> <p>2 PRFM injections administered with 27-G needle: one into the interior vaginal wall, between vagina and urethra, and the other one into the clitoris</p>	<p>Clinical findings from FSDS-R and FSFI patient's questionnaires suggested a trend of improvement after injection. 71% of patients improved from "distressed" to "not distressed", with positive changes in sexual-related dysfunctions</p>	<p>No AE reported.</p> <p>Side effects were observed in 2 patients and including sexual arousal with urination and continuous sexual arousal, ejaculatory orgasm, and spontaneous orgasm. Except for ejaculation, these responses only lasted 1 to 2 weeks. All side effects solved spontaneously.</p>